

PINPOINT LTD

Protecting People at Work



METAL FASCIA CALL BUTTON WITH RESET & DISARM

Backbox: 2 gang 47mm deep (MK877ZIC)

Model No. 132.6

There is a range of Metal Fascia Patient Call Button Units available, which can be configured to different requirements. The standard Unit has an orange Call button and a reassurance LED. To this can be added in any combination, a reset Key Switch, a disarm Key Switch, a combined reset and disarm Key Switch, a ceiling Pull Cord or a 'wander lead' socket. Call Button Units are configured either as 'slaves' to an Advanced Receiver Unit or as separate addressable units. They fit a double-gang back-box and have a brushed-steel effect metal fascia.

An single-gang Metal Fascia Call Button Unit is also available (Model 133.1). It has a call button and a reassurance LED and is installed by connecting to an Advanced Receiver Unit.

Pinpoint's Patient Call options for the 5000 Advanced System do not constitute a Bedhead Service or a Nurse Call System as defined in HTM2015 or NHS Model Specification C49.